


ZEISS FILTER SET 112 HE LED – PENTA-BAND

Use LED 385 (385/30nm) – for Alexa Fluor 350, Alexa Fluor 405, DAPI, Hoechst, etc.









Use LED 475 (469/38nm) - for Alexa Fluor 488, FITC, Cy2, DyLight 488, etc.

Filter Set 112 HE LED

Item Number 489112-9110-000








Use LED 555 (555/30nm) – for Alexa Fluor 546, 568, 594, TRITC, Texas Red, Cy3, etc.


Filter Set 112 HE LED

Item Number 489112-9110-000






Additional components can affect the efficiency. Please also note that the illustration is average.



Use LED 630 (631/33nm) for Alexa Fluor 647, Alexa Fluor 633, Cy5, etc.


Filter Set 112 HE LED

Item Number 489112-9110-000





Additional components can affect the efficiency. Please also note that the illustration is average.




Use LED 735 (735/40nm) – for Alexa Fluor 647, Alexa Fluor 633, Cy5, etc.

Filter Set 112 HE LED

Item Number 489112-9110-000

1 Important Note

This filter set can only be used with Colibr 5 or 7. It does not contain the listed excitation filters. Those are part of the respective Colibri light source.



1 Important Note

Additional components can affect the efficiency. Please also note that the illustration is average.

